


## Catalog Number

95025-E-PE06A

**Vials, Headspace**, Crimp Top, Glass. Amber. 6ml, 20mm Serum Finish, 22x38mm outer dimensions, Beveled Top and Rounded Bottom with a Tapered Neck. Perkin Elmer Compatible. BASIK Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized and considered low actinic. Use with aluminum or steel crimp caps for proper sealing.

Compatible with Perkin Elmer®, Tekmar® and other Headspace autosamplers.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$66.95**

Printed On 2024/05/08