

Catalog Number

95025-E-PE06C

Vials, Headspace, Crimp Top, Glass. Clear. 6ml, 20mm Serum Finish, 22x38mm outer dimensions, Beveled Top and Rounded Bottom with a Tapered Neck. Perkin Elmer Compatible. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass. Use with aluminum or steel crimp caps for proper sealing.

Compatible with Perkin Elmer®, Tekmar® and other Headspace autosamplers.

These vials are also used with Karl Fisher titrators by Mettler and Metrohm.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$59.36

Printed On 2024/05/07