


## Catalog Number

95025-SB

**Storage Box for Headspace Vials.** Blue, polypropylene holds 25 vials. Complete with lid. MICROSOLV Brand. 1 EA.



### Additional Info

These stackable, Storage Boxes are often used for 20mm Headspace Vials for storage or transportation. The Alphanumeric Coding on the boxes makes it easier for Sample Identification. The "egg crate" separators make it easy to keep vials apart.

Supplied with a translucent plastic Lid that keeps condensation from forming on the Septa and minimizes dust or dirt..

- They are Chemically Resistant and can be Autoclaved.
- Temperature Range: -80° to 100°C

**Price**  
**\$70.50**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$70.50**