

## Catalog Number

9502C-0-20CP

**Vial & Cap Kit** includes 100 EA 2ml, Crimp Top, Clear Autosampler Vials and 100 EA matching Crimp Caps with Red PTFE / White Silicone Rubber / Red PTFE Septa, fitted in the Silver Caps. BASIK Brand Easy



### Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the Crimp Top Vials and Caps:

Supplied with Vials: Cat. No [9502C-0-0CV](#) and Caps: Cat No. [9502C-0-20-S](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$84.54**

Printed On 2024/04/29