

## Catalog Number

9502C-0-S

**Vials, Crimp Top, Glass.** Clear, 2ml. An 11mm Serum Finish with 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. BASIK Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, Borosilicate Glass.

Use with 11mm aluminum crimp caps. These narrow opening vials have serum finish tops that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.

Chosen for volatile compounds in GC.

*Lowest cost in class for budget minded laboratories.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$30.31**

Printed On 2024/05/06