

Catalog Number

9502C-0-WCV

Vials, Crimp Top, Glass. Clear, 2ml, with Write-On Patch and fill lines. An 11mm Serum Finish with 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A, pharmaceutical grade, Borosilicate Glass. These wide mouth, vials have a Write-On Patch and fill lines that are epoxy paint for reduced smearing with solvents.

This is our Most Popular Crimp Top Vial with a serum finish.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$35.34

Printed On 2024/05/07