

## Catalog Number

9502C-10CN-M

**Caps, Snap Top**, with Silicone Rubber / PTFE Septa, fitted in smooth, polyethylene, Natural Translucent Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv Brand. 1000/CS.



### Additional Info

*This cap is our original Snap Cap and has been on the market for 30 years. Click [HERE](#) for smaller package quantities..*

*Natural color caps.*

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$315.98**

Printed On 2024/05/02