


Catalog Number

9502C-10CY-M

Caps, Snap Top, with Silicone Rubber / PTFE Septa, fitted in smooth, polyethylene, Yellow Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv Brand. 1000/CS.



Additional Info

This cap is our original Snap Cap and has been on the market for 30 years. Case quantities are the most economical value.

Click [HERE](#) for smaller package quantities.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$341.34

Printed On 2024/05/07