


Catalog Number

9502C-11-MB

Caps, Snap Top, with White Silicone Rubber / Tan PTFE Septa, fitted in smooth, polyethylene, Blue Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. AQ Brand. 100/PK.



MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$45.49

Printed On 2024/04/29