

## Catalog Number

9502C-30CB-M

**Caps, Snap Top**, with Silicone Rubber / PTFE Pre-Slit Septa, fitted in smooth, polyethylene, Blue Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv Brand. 1000/CS.



### Additional Info

This is one of the most popular snap caps on the market and has been available for over 25 years. Preslit septa are normally recommended for Waters Alliance instruments. *Buying in case quantities can save money and time.*

Click [HERE](#) for smaller package quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$427.63**

Printed On 2024/04/28