


Catalog Number

9502C-30CB-X

Caps, Snap Top, with Silicone Rubber / PTFE Pre-Slit Septa, fitted in smooth, polyethylene, Blue Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv Brand. 10,000/CS.



Additional Info

Click [HERE](#) for case quantities of 1,000/CS

Click [HERE](#) for smaller package quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$4,274.40

Printed On 2024/05/05