


Catalog Number

9502C-30CY-M

Caps, Snap Top, with Silicone Rubber / PTFE Pre-Slit Septa, fitted in smooth, polyethylene, Yellow Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv Brand. 1000/CS.



Additional Info

This is one of the most popular snap caps on the market and has been available for over 25 years. *Pre-Slit septa are normally recommended for Waters Alliance instruments.*

Click [HERE](#) for smaller package quantities.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$380.36

Printed On 2024/08/06

Additional Info

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$380.36

Printed On 2024/08/06