

## Catalog Number

9502C-31-MB

**Caps, Snap Top**, with White Silicone Rubber / Tan PTFE Pre-Slit Septa, fitted in smooth, polyethylene, Blue Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. AQ Brand. 100/PK.



### Additional Info

This is one of the most popular snap caps on the market and has been available for over 25 years. *Pre-Slit septa are normally recommended for Waters Alliance instruments.*

Compares to Waters® HPLC Certified Snap Caps.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$70.64**

Printed On 2024/05/05