


Catalog Number

9502C-40T-B

Caps, Crimp Top, with Red Rubber / PTFE Septa, fitted in aluminum, Blue Caps. For use with 11mm Crimp or Snap Top, 12x32mm autosampler vials. BASIK Brand. 100/PK.



Additional Info

The Rubber is softer than Natural Red Rubber and provides better sealing after injection. Great for GC.

Click [HERE](#) for case quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$21.70

Printed On 2024/04/29