


Catalog Number

9502C-41S

Caps, Crimp Top, with 0.05" thick, Silicone Rubber / PTFE, ultra low bleed Septa, fitted in aluminum, Silver Caps. For use with 11mm Crimp or Snap Top, 12x32mm autosampler vials. BASIK Brand. 100/PK.



Additional Info

Ultra low bleed Septa is a popular choice for GC and when siloxanes from silicone rubber may be a problem but you still want to use PTFE. Also used in HPLC when a Crimp Cap is required.

Click [HERE](#) for case quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$27.39

Printed On 2024/05/07