


## Catalog Number

9502C-41S-M

**Caps, Crimp Top**, with 0.05" thick, Silicone Rubber / PTFE, ultra low bleed Septa, fitted in aluminum, Silver Caps. For use with 11mm Crimp or Snap Top, 12x32mm autosampler vials. BASIK Brand. 1000/CS.



### Additional Info

Ultra low bleed Septa is a popular choice for GC and when siloxanes from silicone rubber may be a problem but you still want to use PTFE. Also used in HPLC when a Crimp Cap is required.

Click [HERE](#) for smaller pack quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$246.78**

Printed On 2024/05/05