


Catalog Number

9502C-43T-S

Caps, Crimp Top, with Red Rubber / PTFE Pre-Slit Septa, fitted in aluminum, Sliver Caps. For use with 11mm Crimp or Snap Top, 12x32mm autosampler vials. BASIK Brand. 100/PK.



Additional Info

This Pre-Slit Septa is sometimes used in Solid Phase Micro-Extraction (SPME) methods and when 2ml vials are being used.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$20.94

Printed On 2024/05/03