


Catalog Number

9502C-50CG

Caps, Snap Top, with White PTFE Septa, fitted in smooth, polyethylene, Green Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv Brand. 100/PK.



Additional Info

This PTFE (only) Septa is for single injections or non-volatile samples as it will not reseal.

Click [HERE](#) for case quantities.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$31.75

Printed On 2024/04/25