


## Catalog Number

9502C-90C-B

**Caps, Push Top**, with thinned penetration point, inserted in smooth, polyethylene, Blue Caps. For use with 11mm Crimp or Snap Top, 12x32mm autosampler vials. MicroSolv Brand. 100/PK.



### Additional Info

This cap is useful for non volatile solvents or when evaporation is not considered a problem.

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$18.70**

Printed On 2024/05/18