

## Catalog Number

9502C-E-15G-M

**Caps, Snap Top**, with 0.04" thick, White Silicone Rubber / Natural PTFE Septa, fitted in smooth, polyethylene, Green Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. BASIK Brand. 1000/CS.



### Additional Info

Natural PTFE is very pure and does not have dyes that could extract from it.

Our most economical Silicone Rubber / PTFE Snap Cap. Click [HERE](#) for smaller pack quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$203.41**

Printed On 2024/04/29