

Catalog Number

9502C-E-15R

Caps, Snap Top, with 0.04" thick, White Silicone Rubber / Natural PTFE Septa, fitted in smooth, polyethylene, Red Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. BASIK Brand. 100/PK.



Additional Info

Natural PTFE is very pure and does not have dyes that could extract from it.

Our most economical Silicone Rubber / PTFE Snap Cap. Click [HERE](#) for case quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$28.59

Printed On 2024/05/04