

Catalog Number

9502C-E-31CB-M

Caps, Snap Top, with 0.04" thick, White Silicone Rubber / Yellow PTFE Pre-Slit (Star Shape) Septa, fitted in smooth, polyethylene, Light Blue Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. BASIK



Additional Info

Ultra-low Bleed Septa provide less extractables than most septa and is recommended for LCMS. The star shaped (*) slit offers a much greater target area for an autosampler needle to penetrate the septa and offers minimal resistance. It may not reseal as well as other Pre-Slit Septa.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$341.35

Printed On 2024/04/29