

Catalog Number

9502C-E-35K

Caps, Snap Top, with 0.04" thick, White Silicone Rubber / Natural PTFE Pre-Slit Septa, fitted in smooth, polyethylene, Black Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. BASIK Brand. 100/PK.



Additional Info

Pre-Slit Septa are normally recommended for Waters Alliance and Acquity instruments.

Natural Septa does not have any color pigment and will not introduce extractables.

Click [HERE](#) for case quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$31.76

Printed On 2024/05/18