

## Catalog Number

9502C-E-40T-S

**Caps, Crimp Top**, with Soft-Guard, White Silicone / Red PTFE Septa, fitted in aluminum, Silver Caps. For use with 11mm Crimp or Snap Top, 12x32mm autosampler vials. Soft-Guard Brand. 100/PK.



### Additional Info

Soft-Guard™ Septa provide easier penetration of autosampler needles and better resealing than other brands, for longer instrument up-time and fewer missed injections.

Excellent for HPLC when Crimp Caps are required.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$18.26**

Printed On 2024/04/29