

Catalog Number

9502C-E-WAV

Vials, Snap Top, Glass. Amber, 2ml, with Write-On Patch and fill lines.
An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK



Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized. Amber vials are often used for light sensitive compounds and considered low actinic. Write-On Patch and fill lines are epoxy paint for reduced smearing with solvents.

Lowest cost in class for budget minded laboratories.

Price
\$24.54

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$24.54