


## Catalog Number

9502C-S40-B

**Caps, Snap Top**, with 1.3mm thick, Silicone Rubber / PTFE, ultra-clean Septa, fitted in smooth, soft polyethylene, Light Blue Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv Brand. 100/PK.



### Additional Info

Soft polyethylene caps are more pliable than other Snap Caps and are easier to place on the vial. Popular choice when using a large number of caps at one time. These caps may not seal as well as other Snap Caps and should be used with sample solvents with low volatility.

Ultra-low Bleed Septa provide less extractables than most septa and is recommended for LCMS.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$44.71**

Printed On 2024/05/03