


Catalog Number

9502C-S60-B

Caps, Snap Top, with 1.0mm thick, Orange Rubber / Transparent TEF Septa, fitted in Smooth, soft polyethylene, Blue Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv Brand. 100/PK.



MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$39.54

Printed On 2024/04/30