

Catalog Number

9502C-T-3AP

**Vial & Cap Kit** includes 100 EA 2ml, Snap Top, Amber Autosampler Vials and 100 EA matching Snap Caps with Silicone Rubber / PTFE Pre-Slit, ultra-pure Septa, fitted in the Light Blue Caps. BASIK Brand Easy



### Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 11mm Vials and Caps:

Supplied with Vials: Cat. No [9502C-E-0AV](#) and Caps: Cat No. [9509C-30CB](#)

**Price**  
**\$73.21**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$73.21**