

Catalog Number

9502C-WCV

Vials, Snap Top, Glass. Clear, 2ml, with Write-On Patch and fill lines. An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A pharmaceutical grade, Borosilicate Glass. Write-On Patch is epoxy paint which resists smearing when using solvents.

This is our most popular standard, borosilicate glass Snap Top Vial. For the Most Popular RSA Snap Top Vials, Click [HERE](#).

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$30.61

Printed On 2024/05/04