

## Catalog Number

9502C-WDAV

**Vials, Snap Top, Glass. Amber, 2ml, with Write-On Patch and fill lines. Deactivated. An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 100/PK**



### Additional Info

First hydrolytic, Class A, Type N51A Deactivated / Silanized, pharmaceutical grade, Borosilicate Glass, amberized. Features a wide mouth opening for easy filling. Write-On Patch and fill lines are epoxy paint for reduced smearing.

Amber vials are often used for light sensitive compounds and considered low actinic.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$78.03**

Printed On 2024/05/06