


## Catalog Number

## 9502S-01-B-M

**Caps, Screw Top**, with an Open Hole. "Knurled", polypropylene, Blue Caps. Septa are sold separately and must be inserted manually. For use with 9-425mm thread, 12x32mm autosampler vials. AQ Brand. 1,000/CS.



### Additional Info

These caps are used for our Soft-Guard Brand Caps.

Click [HERE](#) for Septa to use in these caps.

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$45.69**

Printed On 2024/02/25