

## Catalog Number

## 9502S-02CP-K

**Inserts, Low Volume** for Vials, Glass. Clear, 300ul fill volume with a wider, conical point and 6x29mm outer dimensions with plastic spring attached. For use in 2ml autosampler vials. AQ Brand. 1000/CS.



### Additional Info

First hydrolytic, Class A, Type 33, Borosilicate Glass. Features a thicker glass wall and a wider, conical tip for excellent performance to value while minimizing time consuming bubble formation at the tip which can occur with some samples. Polypropylene springs are attached and residual volume is 2.5ul with most autosamplers. This product is manufactured in an ISO 9001:2015 Certified facility.

Use with 9mm Screw Top and 11mm Snap Top 2ml, 12x32mm autosampler vials. Compares to the Polyspring™ Insert.

*This is one of our "Most Popular" Borosilicate Glass Insert and most economical in class. Click [HERE](#) For smaller pack quantities.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$921.49**

Printed On 2024/05/18