

Catalog Number

9502S-02E

Inserts, Low Volume for Vials, Glass. Clear, 250ul fill volume with a conical point and 6x29mm outer dimensions with plastic spring attached. For use in 2ml autosampler vials. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass. 250ul total fill volume and a 100ul usable volume. Features a flat, conical tip for excellent performance. Polypropylene springs are attached and residual volume is 2.5ul with most autosamplers.

Use with 9mm Screw Top and 11mm Snap Top 2ml, 12x32mm autosampler vials.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$114.20

Printed On 2024/05/18