

## Catalog Number

## 9502S-02L-K

**Inserts, Low Volume** for Vials, Glass. Clear, 350ul fill volume with a conical point and 6x31mm outer dimensions. For use in 2ml autosampler vials. MicroSolv Brand. 1000/CS.



### Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass. Features a conical tip for excellent performance. Residual volume is 4ul with most autosamplers.

Spring may be required, sold separately, Can be used in 9mm Screw Top or 11mm Snap Top, 2ml 12x32mm autosampler vials.

*Lowest cost in class for budget minded labs.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$480.05**

Printed On 2024/05/24