

## Catalog Number

9502S-05-B-M

**Caps, Screw Top**, with Silicone Rubber / PTFE Septa, fitted in smooth sided, Robo Style, soft polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. BASIK Brand. 1000/CS.



### Additional Info

Compare to National / Thermo DP™ Screw Caps. Click [HERE](#) for smaller pack quantities.

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$320.98**

Printed On 2024/06/20