

Catalog Number

9502S-05-K-M

Caps, Screw Top, with Silicone Rubber / PTFE Septa, fitted in smooth sided, Robo Style, soft polypropylene, Black Caps. For use with 9-425mm thread, 12x32mm autosampler vials. BASIK Brand. 1000/CS.



Additional Info

Compare to National / Thermo DP™ Screw Caps.

Click [HERE](#) for smaller pack quantities.

Price
\$320.21

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$320.21