

## Catalog Number

## 9502S-05N-K

**Inserts, Low Volume** for Vials, Glass. Clear, 350ul fill volume with a flat bottom and 6x31mm outer dimensions. For use in 2ml autosampler vials. AQ Brand. 1000/CS.



### Additional Info

First hydrolytic, Class A, Type 33, Borosilicate Glass. Features a flat bottom. Residual volume is 6ul with most autosamplers. This product is manufactured in an ISO 9001:2015 Certified facility.

Use with 9mm Screw Top or 11mm Snap Top 2ml, 12x32mm autosampler vials. Lower cost for budget minded laboratories.

Click [HERE](#) for 100/PK or [HERE](#) for 500/PK quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$245.96**

Printed On 2024/05/18