


Catalog Number

9502S-05ND

Inserts, Low Volume for Vials, Glass. Clear, Deactivated, 350ul fill volume with a flat bottom and 6x31mm outer dimensions. For use in 2ml autosampler vials. AQ Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type 33, Deactivated / Silanized Borosilicate Glass. Features a flat bottom. Residual volume is 6ul with most autosamplers. This product is manufactured in an ISO 9001:2015 Certified facility.

Use with 9mm Screw Top or 11mm Snap Top 2ml, 12x32mm autosampler vials. Lower cost for budget minded laboratories.

Click [HERE](#) for case quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$86.75

Printed On 2024/07/27