

Catalog Number

9502S-0AV-K

Vials, Screw Top, Glass. Amber, 2ml, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. BASIK Brand. (10 packs of 100/PK). 1,000/CS.



Additional Info

First hydrolytic, Class A, pharmaceutical grade, Borosilicate Glass, amberized.

Features a wide mouth opening for easy filling and ND9 threads. Lowest cost in class.

Case quantity is the most economical value.

Price

\$317.61

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$317.61