

Catalog Number

9502S-0DCV

Vials, Screw Top, Glass. Clear, 2ml, Deactivated, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A, Deactivated / Silanized pharmaceutical grade, Borosilicate Glass.

Features a wide mouth opening for easy filling. ND9 threads. Lowest cost in class.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$61.18

Printed On 2024/06/20