


## Catalog Number

9502S-10C-R

**Caps, Screw Top**, with Silicone Rubber / PTFE, low durometer Septa, bonded in "knurled", polypropylene, Red Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK.



### Additional Info

The Septa is bonded to the cap for increased performance of penetration that prevents "septa plunge" and prevents them from falling out or dislodging in the cap during transportation..

Very popular and useful for HPLC or GC; save money compared to "Certified" brands.

Click [HERE](#) for case quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$48.44**

Printed On 2024/05/05