


Catalog Number

9502S-10M-G

Caps, Screw Top, with White Silicone Rubber / White PTFE Septa, fitted in "knurled", polypropylene, Green Caps. Recommended for Mass Spec. For use with 9-425mm thread, 12x32mm autosampler vials. Compares to



Additional Info

Made with ultra-pure silicone/PTFE septa and recommended for LCMS or GCMS.

Price
\$54.72

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$54.72