

Catalog Number

9502S-11-MB

Caps, Screw Top, with Tan Silicone Rubber / White PTFE Septa, bonded in "knurled", polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. AQ Brand. 100/PK.



Additional Info

Ultra-pure Septa and is a good choice for LCMS, GCMS & LCMS/MS. this cap is also recommended for Waters® Acquity™ autosamplers.

Click [HERE](#) for case quantities.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$46.46

Printed On 2024/02/25