


## Catalog Number

## 9502S-12Y-M

**Caps, Screw Top, with Silicone Rubber / PTFE, low durometer Septa, bonded in "knurled", polypropylene, Yellow Caps. For use with 9-425 thread, 12x32mm autosampler vials. MicroSolv Brand. 1000/CS.**



### Additional Info

Low durometer makes the Septa easier to penetrate by the autosampler needles.

Click [HERE](#) for smaller pack quantities.

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$468.65**

Printed On 2024/02/25