

## Catalog Number

9502S-30C-R-M

**Caps, Screw Top**, with Silicone Rubber / PTFE Pre-Slit Septa, bonded in "knurled", polypropylene, Red Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand. 1000/CS



### Additional Info

Septa bonded to cap for increased performance of penetration and to prevent "Septa Plunge" into vials during injection. Designed for Waters® Alliance & Acquity instruments.

Click [HERE](#) for smaller pack quantities.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$496.10**

Printed On 2024/05/01