


Catalog Number

9502S-30M-B

Caps, Screw Top, with White Silicone Rubber / White PTFE Pre-Slit Septa, fitted in "knurled", polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK



Additional Info

Pre-slit for easy penetration by autosampler needles. Recommended for Waters® Alliance & Acquity instruments.

Ultra-pure Silicone Rubber and "Cast" PTFE for excellent results in LCMS/MS compares to LCMS Certified brands.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$77.31

Printed On 2024/02/25