

## Catalog Number

## 9502S-31-MB-K

**Caps, Screw Top, with Tan Silicone Rubber / PTFE Pre-Slit Septa, bonded in "knurled", polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials, AQ Brand. 1000/CS.**



### Additional Info

Pre-slit for easy penetration by autosampler needles and to prevent "Septa Plunge". Recommended for Waters® Alliance & Acquity instruments. Ultra-pure Silicone Rubber is a natural Tan color and will not extract pigments. Use this Cap for excellent results in LCMS/MS. Compare to "HPLC Certified" brands.

Click [HERE](#) for smaller pack quantities.

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$629.54**

Printed On 2024/04/25