

## Catalog Number

9502S-33C-B

**Caps, Screw Top**, with Blue Silicone Rubber / White PTFE Pre-Slit Septa, fitted in "knurled", polypropylene, Blue Caps with 6mm holes. For use with 9-425 thread, 12x32mm autosampler vials. MicroSolv Brand.



### Additional Info

Larger, 6mm hole in the cap provides a larger target area for autosamplers that might have difficulty and may strike the cap edge of smaller hole caps, during injection.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$74.40**

Printed On 2024/05/07