

Catalog Number

9502S-35-K

Caps, Screw Top, with Silicone Rubber / PTFE Pre-Slit Septa, fitted in smooth sided, Robo Style, soft polypropylene, Black Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK.



Additional Info

Pre-slit for easy penetration by autosampler needles and to prevent "Septa Plunge" into vials. Compares to Thermo / National Scientific DP® Caps.

Click [HERE](#) for case quantities.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$39.58

Printed On 2024/02/25