

## Catalog Number

## 9502S-80C-B

**Caps, Screw Top**, with Natural Rubber / TFE Septa, fitted in "knurled", polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK.



### Additional Info

Natural Rubber does not contain siloxanes like Silicone Rubber. TFE is similar to PTFE when comparing inertness.

Click [HERE](#) for case quantities.

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$38.54**

Printed On 2024/02/25